



# SUPRA G

PORCINE AORTIC VALVE CONDUIT

# SUPRA-G

## The Biological option for Bentall procedure

- Unique porcine aortic valved conduit
- Ultra Low Porosity Woven Polyester
- Offers excellent haemostasis



- Preassembled device potentially reduces:
  - Bypass time
  - Surgeon labor
  - Procedure time
  - Cross-clamp time


# Product Information

## Labcor TLPB-A SUPRA-G

Catalogue Number	Valve Information	Graft Information	
	Valve Size (mm)	Internal Diameter (mm)	Length (mm)
TLPB-A SUPRA-G 19	19S	22	100
TLPB-A SUPRA-G 21	21S	24	
TLPB-A SUPRA-G 23	23S	26	
TLPB-A SUPRA-G 25	25S	28	
TLPB-A SUPRA-G 27	27S	30	
TLPB-A SUPRA-G 29	29S	32	

### Presentation:

TLPB-A Supra G is presented sterile and non-pyrogenic. Immersed in 0.2% glutaraldehyde solution.





**CE**  
0434

**Labcor Laboratórios Ltda**  
Av. Silva Lobo, 1670 - Nova Granada  
Belo Horizonte/MG - Brasil  
CEP 30431-262  
Tel.: 55 31 3371-1371  
Depto Comercial.: 55 21 3722-1600  
contato@labcor.com.br - www.labcor.com.br

**P/N: 99B0091/A**